DATA BASE FOR MONITORING BIOMEDICAL DEVICES

Data Base for monitoring biomedical devices

The aim of our project is to create a data base for the Hospital "La Paz" in order to organize all the information related with the doctors, their patients, and the biomedical devices in use.

The six entities are:

- 1. Appointment
- 2. Devices
- 3. Doctors
- 4. Manufacturer
- Patient

Each entity has different properties, as well as an Id which identifies it. We can observe this in the ER diagram.

USER STORIES:

- As a patient, I need a device to improve my quality of life.
- As a doctor, I treat a patient to monitor their diagnosis.
- As a doctor, I accept an appointment to check on a patient.
- As a patient, I order an appointment with the doctor.
- As a manufacturer, I check what devices do the hospital need in order to fabricate them.

USE CASES:

Order devices:

Actor: Doctor

- ~ Description: The doctor orders an amount of devices in a specific date.
- ~ Preconditions: An appointment must had been taken place between the doctor and the patient. The doctor must had decided that the patient needs a device according to its diagnosis.

Standard path:

- 1. Select doctor
- Select the option order device
- 3. Select the device requested and the amount

• Confirms an appointment :

Actor: Doctor

- ~ Description: The doctor enters the application and requests an appointment on a specific date.
 - ~ Preconditions: The doctor must be registered in the application.

Standard path:

- Select doctor
- 2. Select the option request appointment

• Request appointment :

Actor: Patient

- ~ Description: The patient enters the application and requests for an appointment on a specific date.
 - ~ Preconditions: The patient must be registered in the application.

- 1. Select patient
- 2. Select the option request appointment
- 3. Select the date of the appointment

Register patient:

Actor: Patient

~ Description: The patient enters on the application and introduces his personal information to register.

~ Preconditions:

Standard path:

- Select patient
- 2. Select the option register patient
- 3. Introduce the personal information

• Cancel an appointment :

Actor: Doctor

~ Description: The doctor enters the application and cancels an appointment previously made by a patient for a specific date.

~ Preconditions: The patient must have asked for an appointment.

Standard path:

- Select doctor
- 2. Select the option cancel appointment

View information about device:

Actor: Patient

~ Description: The patient enters the application and checks the specific information about his device.

~ Preconditions: The patient needs to be registered and must have a device.

- 1. Select patient
- 2. Select the option view information about device

View devices needed:

Actor: Manufacturer

- ~ Description: The manufacturer enters the application and edits his personal information.
 - ~ Preconditions: The Hospital must have requested devices.

Standard path:

- 1. Select manufacturer
- Select the option view devices needed

· Modify data:

Actor: Manufacturer

- ~ Description: The manufacturer enters the application and modifies his own information.
- ~ Preconditions: The manufacturer must be registered in the application.

Standard path:

- 1. Select manufacturer
- 2. Select the option edit personal information
- 3. Select the information that he wants to change

• Cancel appointment :

Actor: Patient

- \sim Description: The patient enters the application and cancels the appointment.
- ~ Preconditions: The patient must have requested for an appointment.

- 1. Select patient
- 2. Select the option cancel appointment

• View information about patient :

Actor: Doctor

- ~ Description: The doctor enters the application and checks the information about the patient.
- ~ Preconditions: The patient needs to be registered and they must have had an appointment.

Standard path:

- Select doctor
- 2. Select the option view information about patient

· Modify data:

Actor: Doctor

- ~ Description: The doctor enters the application and modifies his own information.
- ~ Preconditions: The doctor must be registered in the application.

Standard path:


- Select doctor
- 2. Select the option edit personal information
- 3. Select the information that he wants to change

Modify data :

Actor: Patient

- ~ Description: The patient enters the application and modifies his own information.
- ~ Preconditions: The patient must be registered in the application.

- 1. Select patient
- 2. Select the option edit personal information
- 3. Select the information that he wants to change



REQUIREMENTS

FUNCTIONAL:

- 1 The system stores info. about the patient in the D.B.
- 2 The system stores info. about the devices in the D.B.
- 3 The doctor can add a device to the system
- 4 The doctor can modify the patient's data
- 5 The doctor has access to the patient's data
- 6 The doctor can add data to the system
- 7 The doctor can print the list of the patients that he treats
- 8 The doctor can assign a patient to him
- 9 The doctor can modify his own data
- 10 The patient can choose a doctor for his appointment
- 11 The patient can obtain his own data
- 12 The patient can obtain his own device information
- 13 The patient can modify his own data
- 14 The access to the app requires login
- 15 The system has three roles of users

NON-FUNCTIONAL:

- i The system should respond to the user interaction in less than two seconds
- ii The system should be able to support multiple users at the same time without reducing performance
- iii In case of system failure, the system can be able to recover within 30 minutes
- iv The system does regular backups of the D.B.

TRACEABILITY MATRIX:

	Requirements .																			
بر دون دون	25	1	2	3	ц	2	6	7	8	9	10	11	12	13	14	15	I	I	耴	<u> 17</u>
	A		/	\checkmark					\checkmark	✓				\checkmark	,		V	\checkmark	\checkmark	/
	B				✓				V	\checkmark				$\sqrt{}$	✓		V	\checkmark	\checkmark	
	<u>_</u>	V				/	V			V						\checkmark	V	✓	/	/
	D	\checkmark											\checkmark	1			V	V	/	/
	E	\checkmark			\checkmark				√	\checkmark				\checkmark			V	V	/	V
	F	\checkmark	✓				√	\checkmark					\checkmark	\checkmark			\checkmark	V	/	
	6		✓	\checkmark									\checkmark	\checkmark			\checkmark	✓	\checkmark	/
	H		1										$\sqrt{}$	$\sqrt{}$			V	ν	$\sqrt{}$	/
	I	\checkmark				\checkmark	\checkmark			\checkmark			\checkmark			\checkmark	$\sqrt{}$	V	/	V
	J	\checkmark	\checkmark							\checkmark	$\sqrt{}$		$\sqrt{}$	\checkmark			\checkmark	V	/	V
	K								\checkmark		$\sqrt{}$	$\sqrt{}$	\checkmark	\checkmark	\checkmark		V	V	/	~
	ļ	\checkmark				\checkmark	$\sqrt{}$						\checkmark		$\sqrt{}$		\checkmark	\checkmark	/	/
		ı	l		ı				I	I	I	•								

Patients

p_id	Name	Surname	Device	B date	Diagnosis
1.	Maria	Gorgojo	Pacemaker	2003-9-22	Arrhythmia
2	Sandra	Cabrera	Prosthetic Limb	2003-7-25	Amputated Leg

Doctor

do_id	Specialty	Name	Surname
1	Cardiologist	José Luis	Luchelli
2	Traumatologist	Roberto	Rodríguez

Device

d_id	Туре	manufacturer_id
1	Pacemaker	1
2	Insulin pump	2

Manufacturer

m_id	Name	Adress	Phone number
1	Lucia	Andes Av.	689054844
2	Ana	Serrano street	629941093

Appointment

a_id	Date	Description	doctor_id	patient_id
1	2024-6-2	Revision	1	1
2	2024-7-5	Revision	2	2

Device_doctor: Orders

doctor_id
FK to do_id
1
2

Device_patient

device_id	patient_id
FK to d_id	FK to p_id
1	1
2	2

